




ZIMMER BIOMET

Your progress. Our promise.®



Spine Solutions

Gallery™ Laminoplasty Spine System

Comprehensive and
intuitively designed

Gallery Laminoplasty

SPINE
SYSTEM

The Gallery Laminoplasty Spine System is a comprehensive and intuitive system designed to prevent spinal cord impingement in the lower cervical and upper thoracic spine.



FLEXIBILITY

- Cobra head design
- Variety of plate styles and sizes
- Self-tapping and self-drilling screws



SIMPLICITY

- Design incorporates quick, secure connection between the screws and inserter
- Spin Tight Technology makes screw insertion easy

SMART PLATE DESIGN

The Gallery Laminoplasty Spine System offers two plate styles. Both cobra head plate designs are low profile and offer multiple points of fixation for ideal screw placement.

- Design incorporates a quick, secure connection between the screws and inserter
- Spin Tight Technology makes screw insertion easy

MULTIPLE IMPLANT OPTIONS

VARIETY OF PLATE SIZES AND STYLES

The Gallery Laminoplasty Spine System plates are available in hook or standard configuration, in a variety of lengths from 8 mm to 16 mm (2-mm increments).



Standard plate

Plate with hook



SELF-TAPPING AND SELF-DRILLING SCREWS

Gallery Laminoplasty Spine System screws are available in self-drilling and self-tapping options in multiple lengths from 5 mm–11 mm (2-mm increments).


Screws are available in 2.4 mm and 2.8 mm diameters.



Self-tapping Screws



Self-drilling Screws



ERGONOMIC INSERTER

QUICK, SECURE CONNECTION

The design behind the Gallery Laminoplasty Spine System inserter ensures a quick, secure interface between the screw and inserter.

- Quick, secure connection
- Wheel mechanism
- Ergonomic handle fits in the palm of your hand

SPIN TIGHT TECHNOLOGY

Spin Tight Technology offers a wheel mechanism on the handle, the location of which allows the surgeon to simultaneously place pressure on the inserter with the palm of the hand and drive the screw with minimal effort.





ZIMMER BIOMET
Your progress. Our promise.®

800.447.3625/zimmerbiomet.com

©2017 Zimmer Biomet Spine, Inc. All rights reserved.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet Spine, Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet Spine. This material is intended for health care professionals, the Zimmer Biomet Spine sales force and authorized representatives. Distribution to any other recipient is prohibited.

For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert and www.zimmerbiomet.com.

0436.1-US-en-REV1017